

# Adjustable Knee Brace

## 可調膝關節支架

**Model: TD-TX**

### Features:

- Adjustable chuck, gradual recovery angle adjustment, simple and convenient, can walk and exercise
- Multiple non-slip straps + high-viscosity Velcro, tightly integrated with the legs, loose and adjustable
- Comfortable non-slip lining, not easy to fall off when worn  
Suitable for: knee joint injuries and rehabilitation exercises
- 可調節卡盤, 漸進式恢復角度調節, 簡單方便, 可行走鍛煉
- 多條防滑綁帶+高粘魔術貼, 與腿部緊密結合, 寬鬆可調
- 舒適防滑內襯, 穿戴不易脫落
- 適用於: 膝關節部位受傷和康復鍛煉

### Specifications

Product Model 產品型號	TD-TX
Product Function 產品功能	rehabilitation training / fitness 康復訓練 健身
Material 材質	aluminum alloy 鋁合金
Mode 模式	可調節 adjustable
Model 型號	均碼 free size
Suitable For Thigh Circumference 適合大腿圍	49-65cm
Suitable Calf Circumference 適合小腿圍	36-52cm
Product Size 產品尺寸	55*29*5cm
Product Weight 產品重量	933g
Package Size 包裝尺寸	59.5*15*15cm
Package Weight 包裝重量	1.30kg



### Instructions

1. 將支具平鋪展開, 根據固定綁帶的位置, 將墊片從上到下的順序, 固定在腿上
2. 調整好位置後, 卡盤位置正對膝關節, 將支架粘貼在腿部的左右兩側正中
3. 將所有固定綁帶收緊, 並粘扣好即可

#### 卡盤角度調節:

1. “Flex”側調整彎曲角度, 密曲角度在0-120度之間調節
2. “Ext”側調整伸展角度, 伸展角度0-90度之間調節
3. 彎曲角度於伸展角度共同作用, 可以實現卡盤的固定, 範圍活動及自由活動

#### ◆ 注意事項:

1. 最好由醫護人員或者看護人幫助固定
2. 建議先固定好支具後, 再調節角度, 這樣不易牽動傷口
3. 固定綁帶時, 注意不要太緊或太鬆, 太緊影響血液循環, 太鬆護具沒有作用