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
Finger Stretching Trainer

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Instruction Manual

(FT2-B)

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# Contents

Product introduction	02
User guide	04
Safety precautions	11
Safety precaution include	12
Troubleshooting	14
Technical parameters	15
After service	16

After service, Thank you for choosing the finger trainer.

Please read this manual carefully before using this product.

Please keep this manual properly for guidance when using the stretching trainer.

The finger training device has two modes: the passive mode and mirror training mode. Our device aims to maintain and expand the range of motion of finger joints, prevent joint contracture and muscle atrophy, greatly promote the reorganization and compensation of the nervous system, and help patients recover as soon as possible.

The recovery of hand dysfunction after stroke is a long-term and arduous process. Patients' conditions, rehabilitation courses, and curative effects are very different. Under normal circumstances, the recovery of hand function is slower than other parts, and long-term rehabilitation training is needed.

If you have any questions in use, please contact: [info@ziglite.com](mailto:info@ziglite.com)

The manual version may upgrade due to changes in software version or technical specifications without prior notice.


# Product Introduction

## Scope Of Application

It is suitable for rehabilitation training in (families) institutions for patients with finger joint recovery after injury or limb dysfunction sequelae caused by brain injury, stroke, cerebral palsy, spinal cord injury, orthopedic surgery, multiple sclerosis, etc.


## Product Composition

### Trainer Host


# Product Introduction


## Training Gloves

## Mirror Glove

## Power Adapter

# User Guide

## Buttons, screen icons in detail

1

### Strength Level Display & Settings

default to 9th gear

Press this button to set the training strength level, and the corresponding function display above will flash on the digital tube. At this time, you can use the "+" and "-" keys to set the range from 1 to 9 gears.

2

### Finger Opening Time Level Display & Settings

default to 1st gear

Press this button to set the finger opening time level, and the corresponding function display above will flash on the digital tube. At this time, you can use the "+" and "-" keys to set the range from 1 to 9 gears.

# User Guide

3

default to 1st gear

## Finger Contraction Time Level Display & Settings

Press this button to set the finger contraction time level, and the corresponding function display above will flash on the digital tube. At this time, you can use the "+" and "-" keys to set the range from 1 to 9 gears.

4

default to 15 minutes

## Training Time Display & Time Setting

Press this button to set the training time. At this time, the corresponding function display, namely the fourth and fifth digital tubes, will flash and can be set using the "+" and "-" keys. The setting range is 5-60 minutes, and the graduation value is 5 minutes.

5

default to training mode

## Training/Mirror Mode Display

The upper figure shows the working status indication of the training mode;

The figure below shows the working status indication of mirror mode.

(All green lights indicate the current mode status)

# User Guide

6

## Return Button

Press this key to return the machine to the default mode of initial startup in any mode.

7

## Mirror Mode Button

Press this button to enter/exit the mirror training mode.

8

## + "Setting Button

Press this button to increase the settings for training strength, finger opening time, and finger contraction time.

9

## - "Setting Button

Press this button to reduce the settings for training strength, finger opening time, finger contraction time, and training time.

10

## Start / Pause button

Press this button to start and pause training settings.

# User Guide

## Use Guide

**1**

Place the therapeutic device host

Place the therapeutic device host on the desktop, not on the wet, soft cloth or unclean desktop.

**2**

Connect air pipe and DC plug

When connecting the training gloves, it is necessary to clip the air pipe connector into the corresponding interface of the host.

**3**

Wear training (mirror) gloves

The guardian shall help the patient wear gloves and pull the glove Velcro to the wrist to bond and tighten.

# User Guide

## Use Guide

4

### Power on

Turn down the "power switch" button of the trainer host to the "ON" position to turn on the power, and the display screen will light up to enter.

5

### Start hand training

Press the "cycle mode / mirror mode" button to start and pause training at any time.

6

### Shut off


After the training, turn the "power switch" button of the trainer host to the "OFF" position to power off.

**Note:** The training mode controls the timing of fingers through program control, which first controls the four fingers outside the thumb to contract first, and the thumb will delay contraction by 2 seconds. The advantage is to avoid causing the thumb and index finger to press against each other when the five fingers contract at the same time, which affects the training effect.

# User Guide

## Glove Wiring Diagram


### Glove Joint Insertion Effect

# User Guide

## Air Valve Switch Control


Each glove trachea of the training gloves is equipped with a single independent air valve switch, which can be turned on or off to carry out finger targeted selection training.

### Open The Air Valve

The direction of the air valve is the same as that of the air pipe. The air valve is open, and the corresponding fingers can undergo rehabilitation training.

The air valve switch is on by default.

### Close The Air Valve

The direction of the air valve and the air pipe is cross shaped, the air valve is closed, and the corresponding fingers are not trained.

# Safety Precautions

**[warning] failure to comply with the following safety instructions may lead to serious accidents such as fire, electric shock and casualties.**

- Do not disassemble, repair, or modify the components of the equipment by yourself to avoid abnormalities in the treatment device.
- Do not store components in or near water heaters, microwave ovens, cooking equipment, or high-pressure containers to avoid fire caused by overheating.
- It cannot be sterilized or autoclaved.
- Do not continue to use it when the alarm prompt appears.
- It is suitable for hand circulation and mirror movement rehabilitation training. Please do not use it for other purposes to avoid accidents.
- Do not store the parts in a place that is too hot or too cold to avoid accidents.
- The validity period of the product is 5 years, and the disposal of the product at the end of its service life should comply with the requirements of local laws and regulations.
- Please keep the power connector dry and water-free, please do not plug or unplug the power supply with wet hands to avoid electric shock or injury.

## **Working Environment**

- 1) Ambient temperature: 10°C~35°C
- 2) Relative humidity: 30%~75%, excluding condensation
- 3) Atmospheric pressure range: 700hpa~1060hpa
- 4) Power supply: AC100~240V, 50/60Hz (main unit power supply of the treatment device DC12V)

## **Storage And Transportation**

- During transportation, please fix the trainer in a relatively stable place to prevent collision and damage to the trainer.
- Please prevent liquids or foreign objects from entering the trainer to avoid damaging internal electronic components and causing malfunctions.
- Please do not place it in places with strong sunlight or ultraviolet radiation for a long time to avoid causing discoloration of the shell.
- Do not place other items on the product when saving it.
- Do not fold the gloves forcefully to avoid damaging the soft-actuated actuator and causing malfunctions.
- Do not store the trainer in a humid environment or strong magnetic field interference that is not suitable for working with electronic devices.

# Safety Precaution Include

- In case of skin allergy and other discomfort during the use of the product, please stop using and consult a doctor.
- For patients with diseases or wounds or being treated for them, or with previous medical history, please consult a doctor after being familiar with the following contents.
- In case of sprain, trauma, bone injury of unknown condition, severe osteoporosis, untreated fracture or people who need bracing after treatment, do not use this product without professional confirmation.

## Contraindication

- Patients with severe blood flow disorder such as thoracic venous tumor and hemocolumn disease;
- Patients with ulcer, wound, purulent fruit, bleeding, scald, inflammation, etc. on the affected hand.

## Other Precautions

- 1) In case of abnormal experience or body feeling during use, stop using it and use it under the guidance of a doctor after seeing a doctor;
- 2) In case of redness, pruritus and rash during use, stop use and seek medical attention;
- 3) Do not use it together with other equipment or in the state of applying drugs;
- 4) Do not use it alone for unconscious people or young children;
- 5) Do not use this product for any purpose other than the purpose of use.

## Before Use

- 1) Before wearing training gloves, please try to run the therapeutic device host for 1 minute to confirm that all parts function normally;
- 2) Patients with mandatory systolic spasm need to massage before using the therapeutic device;
- 3) Note: before training, please confirm that the non muscle tension of the hand to be trained is too high (above grade II).

# Safety Precaution Include

## **In Use**

In case of the following conditions during operation, please stop using immediately and disconnect the power supply.

- 1) The trainer host emits abnormal sound;
- 2) The surface temperature of the trainer host is too high;
- 3) Water ingress into the trainer host;
- 4) Users experience pain, swelling, or other discomfort symptoms;
- 5) If any spasms occur during the training process, immediately stop using it.

## **Usage Time**

- 1) Suggest selecting appropriate training methods based on one's actual situation;
- 2) The duration of each use should not be too long. If there is soreness in the hands the next day after training, it is necessary to reduce the number and duration of daily training;
- 3) It is recommended to use it 2-4 times a day, each time for about 15 minutes, with a minimum interval of 45 minutes between uses;
- 4) If it is necessary to replace or extend the use time during use, please consult professional medical personnel;
- 5) Long term continuous use can easily lead to high internal temperature of the machine, which may cause phenomena such as sleep or asynchronous movement. After the machine is fully cooled, it will automatically return to normal.

## **Use Of Power Adapter**

- 1) The power adapter must be plugged and unplugged by a patient's supervisor;
- 2) After using the trainer, please disconnect the power adapter from the 220V power supply;
- 3) Do not use a power adapter that is not equipped with this machine.

# Troubleshooting

Problem	The Reason	Measures
Flip the "power switch" down to turn on the machine	Poor power connection	Check the power cord and power adapter, the connection between the power adapter and the host
Training gloves are not working	Training gloves not properly connected to the host or damaged	Check the physical status of the training gloves and the connection between the host and the training gloves
The treatment gloves are not installed or removed from the treatment device	Training glove tracheal joint not rotated or twisted out	Rotate the training glove tracheal joint counterclockwise and unscrew it

**Note:** If you still can't use the above method, please contact us.

## Care And Maintenance

### Clean And Disinfect

- Do not wash training gloves and mirrored gloves directly with water. Dip a soft cloth in water or neutral detergent, wring it out, and wipe the gloves.
- When dust or dirt appears on the surface of the trainer host, users can use medical alcohol with a concentration of 75% to wipe it. When wiping, please dip a small amount of alcohol into a dry cloth to avoid allowing alcohol to drip or flow into the product.
- Clean with clean and dry soft cloth, and dry in a cool and ventilated place.

### Note

1. The use of cleaning and disinfecting agents other than alcohol may damage the product, shorten the service life or cause safety risks.
2. It is forbidden to use high-temperature drying and exposure to the sun for quick drying or disinfection.
3. Do not use gasoline, thinner or similar organic solvents to clean any part of this product. Do not use radiation, steam, ethylene oxide and other methods for disinfection, otherwise the product will be damaged.

# Technical Parameters

## Product Specifications


Product Name	Finger Stretching Trainer	Waterproof Level	IPX0
Product Model	FT2-B	Velocity Setting	1-9
Security Type	I Type	Time Setting	5-60 minutes
Power Supply	AC100-240V 50/60Hz	Running Noise	≤60dB
Host Size	251*204*111mm	Training Mode	Cycle training/ Mirror training
Power	12W	Standby	1W

## Glove Size

Size	S	M	L	XL	XXL
Hand length (mm)	145-165	160-185	180-195	195-210	210-235
Hand width (mm)	65-80	80-95	90-105	100-115	115-125
Middle finger circumference (mm)	40-50	50-60	60-70	70-90	90-100

\* The corresponding model can be selected according to the actual hand measurement of the user. If two or more measured values are at the critical point, it is recommended to select one size larger.

\* Women generally choose s/m, men generally choose l/xl, and the larger hand is XXL.

# After service

After sales email : [info@ziglite.com](mailto:info@ziglite.com)

After sales address : Room B3, 18/F Bonsun Industrial Building,  
364-366 Sha Tsui Road, Tsuen Wan, Hong Kong

Packing List		
Number	Name	Quantity
1	trainer host	1 set
2	training gloves	1 set
3	wireless mirror gloves	1 set
4	power adapter power cord	1 set
5	product manual	1 set
6	USB data cable	1 set

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Office Address : Room B3, 18/F Bonsun Industrial Building, 364-366 Sha Tsui Road, Tsuen  
Wan, Hong Kong

Email : [info@ziglite.com](mailto:info@ziglite.com)

Fax: +852 3544 0462

Office : +852 3107 7500

Whatsapp : +852 6827 2010

NO. FT2-B (CH-EN)

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Ziglite Smart Health Care Products Co., Ltd.