

Aluminum alloy walking aid with pulley

帶輪鋁合金助步器-彎曲扶手

Model: LR-RP

Features:

- The product is easy to operate, foldable, convenient to travel and carry, and the height can be adjusted at will
- Aluminum alloy material, safe and stable, the handle is wrapped with sponge, so it won't hurt your hands
- Upgrading and bolding the H crossbar to carry a higher load-bearing capacity of up to 100 kg
- S-bend design for assisting to stand up, no need for help from others
- The design of rear caster wheels makes walking easier
- Applicable people: the elderly, those with inconvenient legs and feet, and those who have rehabilitated legs and feet

- 產品簡易好操作,可折疊,方便出行攜帶,高度可隨意調節
- 鋁合金材質,安全穩固,手把用海綿包裹,不傷手
- 升級加粗H橫杆承重更高承重可達100公斤
- 輔助起身S彎設計,自起身不用別人攙扶
- 後腳滾輪設計,助行更方便
- 適用人群:老年人,腿腳不方便者,腿腳復健者

Specifications

Product Model 產品型號	LR-RP
Product Material 產品材質	Aluminum Alloy 鋁合金
Expanded size 展開尺寸	77*51*43cm
Folding size 折疊尺寸	77*51*10cm
Adjustable height 可調節高度	77-87cm
Packing size 包裝尺寸	80*56*9cm
Package weight 包裝重量	2.5KG



Instructions

1. Take out the product and open both sides, press and hold the button in the middle of the top of the product when closing, the four feet can be adjusted at will
2. Remove the front telescopic rod and hold down the button, and the multi-alignment hole position of the lenticular rod is stuck in and the replacement is completed.

- 1.將產品取出打開兩邊即可,收的時候按住產品頂部中間的那個按鈕即可,四腳可隨意調節高低
- 2.將前面的伸縮桿按住按鈕拆下,倫桿多準孔位卡進去即更換完成